


# Mono-2-O-(p-toluenesulfonyl)-β-cyclodextrin

货号: B27483



产品信息

|      |                                                                                                                                                                       |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 生物活性 | Mono-2-O-(p-toluenesulfonyl)-β-cyclodextrin is a biochemical reagent that can be used as a biological material or organic compound for life science related research. |
| CAS  | 84216-71-7                                                                                                                                                            |
| 中文名称 | 单-2-O-(对甲苯磺酰)- β-环糊精水合物                                                                                                                                               |
| 分子量  | 1289.17                                                                                                                                                               |
| 体外研究 |                                                                                                                                                                       |
| 体内研究 |                                                                                                                                                                       |
| 形式   |                                                                                                                                                                       |
| 运输条件 | Room temperature in continental US; may vary elsewhere.                                                                                                               |
| 保存条件 | Please store the product under the recommended conditions in the Certificate of Analysis.                                                                             |
| 溶解性  |                                                                                                                                                                       |
| 纯度   |                                                                                                                                                                       |