

Ganglioside sialidase (AuSialidase M2)

货号: B27484

Ganglioside sialidase(AuSialidase M2)

产品信息

生物活性	Ganglioside sialidase (AuSialidase M2) from <i>Arthrobacter ureafaciens</i> . Ganglioside sialidase is a highly specific N-acetylneuraminidase. Ganglioside sialidase can hydrolyze the internal sialic acid of GM1 under optimal condition with sodium cholate.
CAS	2170331-76-5
中文名称	
分子量	
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	