

# Adenylate

货号: B27532

## Adenylate

产品信息

生物活性	Adenylate (Adk), i.e., myokinase, is a phosphotransferase enzyme, is often used in biochemical studies. A denylate catalyzes the interconversion of adenosine phosphates. Adenylate monitors phosphate nucleoti de levels inside the cell, it plays an important role in cellular energy homeostasis.
CAS	9013-02-9
中文名称	
分子量	
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	