


Guanosine 5'-triphosphate-5'-adenosine

货号: B27547



产品信息

生物活性	Guanosine 5'-triphosphate-5'-adenosine, the 5' cap analog , is a fluorescent substrate analog.
CAS	10527-47-6
中文名称	
分子量	772.41
体外研究	<p>Guanosine 5'-triphosphate-5'-adenosine (GpppA) labeled with pyrene at the 3'-O position of adenosine acts as an artificial substrate.</p> <p>Fluorescently labeled Guanosine 5'-triphosphate-5'-adenosine (GpppA) and GpppG analogs as potential substrates that represent a reasonable compromise between the structural complexity and requirements of the enzyme.</p> <p>The accuracy of these methods have not been independently confirmed. They are for reference only.</p>
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	