


# Heparin disaccharide I-S tetrasodium

货号: B27738



## 产品信息

|      |                                                                                                               |
|------|---------------------------------------------------------------------------------------------------------------|
| 生物活性 | Heparin disaccharide I-S tetrasodium is the major disaccharide produced from heparin by heparanases I and II. |
| CAS  | 136098-10-7                                                                                                   |
| 中文名称 | 肝素二糖I-S,四钠盐                                                                                                   |
| 分子量  | 665.40                                                                                                        |
| 体外研究 |                                                                                                               |
| 体内研究 |                                                                                                               |
| 形式   |                                                                                                               |
| 运输条件 | Room temperature in continental US; may vary elsewhere.                                                       |
| 保存条件 | Please store the product under the recommended conditions in the Certificate of Analysis.                     |
| 溶解性  |                                                                                                               |
| 纯度   |                                                                                                               |