


# Brilliant green

货号: **B27771**

## 产品信息

生物活性	Brilliant green is a cationic dye used to color silk and wool. Brilliant green inhibits propagation of mold, intestinal <b>parasites</b> and <b>fungus</b> . Brilliant green is effective against Gram-positive <b>bacteria</b> .
CAS	633-03-4
中文名称	钻石绿
分子量	482.63
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	