

Phospholipase A2

货号: B27802

Phospholipase A₂

产品信息

生物活性	Phospholipase A2 (PLA2) catalyzes the hydrolysis of the sn-2 position of membrane glycerophospholipid s to liberate arachidonic acid (AA). Phospholipase A2 is a member of the class of heat-stable, calcium-d ependent enzymes, is often used in biochemical studies.
CAS	9001-84-7
中文名称	
分子量	
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	