


# Lithium 3,5-diiodosalicylate

货号: **B27807**

## 产品信息

|      |                                                                                                                                                        |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------|
| 生物活性 | Lithium 3,5-diiodosalicylate is a biochemical reagent that can be used as a biological material or organic compound for life science related research. |
| CAS  | 653-14-5                                                                                                                                               |
| 中文名称 | 2-羟基-3,5-二碘苯甲酸 锂盐                                                                                                                                      |
| 分子量  | 395.85                                                                                                                                                 |
| 体外研究 |                                                                                                                                                        |
| 体内研究 |                                                                                                                                                        |
| 形式   |                                                                                                                                                        |
| 运输条件 | Room temperature in continental US; may vary elsewhere.                                                                                                |
| 保存条件 | Please store the product under the recommended conditions in the Certificate of Analysis.                                                              |
| 溶解性  |                                                                                                                                                        |
| 纯度   |                                                                                                                                                        |