


货号: B28047



| | |
|------|---|
| 生物活性 | DODAP (hydrochloride) is an ionizable lipid. DODAP (hydrochloride) has the potential for the research of gene delivery. |
| CAS | 329009-00-9 |
| 中文名称 | |
| 分子量 | 684.52 |
| 体外研究 | |
| 体内研究 | |
| 形式 | |
| 运输条件 | Room temperature in continental US; may vary elsewhere. |
| 保存条件 | Please store the product under the recommended conditions in the Certificate of Analysis. |
| 溶解性 | |
| 纯度 | |