


# ZnAF-2

货号: **B28109**

## 产品信息

生物活性	ZnAF-2, 6-iso- is a cell-impermeable compound that acts as a high-affinity Zn <sup>2+</sup> -specific fluorescent probe (K <sub>d</sub> = 2.7 nM). It exhibits low basal fluorescence with an approximately 51-fold increase in fluorescence intensity upon stoichiometric (1:1) binding to Zn <sup>2+</sup> . Little affinity for Ca <sup>2+</sup> , Mg <sup>2+</sup> , Na <sup>+</sup> or K <sup>+</sup> .
CAS	321859-11-4
中文名称	
分子量	572.61
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	