

Mid-chain triglyceride of caprylic/caprolic acid

货号: B28228

Mid-chain triglyceride of caprylic/caprolic acid

产品信息

| | |
|------|---|
| 生物活性 | Mid-chain triglyceride of caprylic/caprolic acid can be used as a co-solvent. |
| CAS | |
| 中文名称 | |
| 分子量 | |
| 体外研究 | |
| 体内研究 | |
| 形式 | Liquid |
| 运输条件 | Room temperature in continental US; may vary elsewhere. |
| 保存条件 | |
| 溶解性 | |
| 纯度 | |