

# Tyrosinase

货号: **B28244**

## Tyrosinase

产品信息

生物活性	<b>Tyrosinase</b> (Polyphenol oxidase) is a rate-limiting enzyme that controls the production of melanin and is encoded by TYR gene. <b>Tyrosinase</b> is mainly found in melanosomes synthesized by skin melanocytes.
CAS	9002-10-2
中文名称	
分子量	
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	