


# Diosgenin-3-O-β-D-sodium glucuronide

货号: **B28334**



产品信息

|      |                                                                                                                                                                |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 生物活性 | Diosgenin-3-O-β-D-sodium glucuronide is a biochemical reagent that can be used as a biological material or organic compound for life science related research. |
| CAS  | 32471-99-1                                                                                                                                                     |
| 中文名称 | 薯蓣皂苷元-3-O-B-D-葡萄糖醛酸钠盐                                                                                                                                          |
| 分子量  | 612.73                                                                                                                                                         |
| 体外研究 |                                                                                                                                                                |
| 体内研究 |                                                                                                                                                                |
| 形式   |                                                                                                                                                                |
| 运输条件 | Room temperature in continental US; may vary elsewhere.                                                                                                        |
| 保存条件 | Please store the product under the recommended conditions in the Certificate of Analysis.                                                                      |
| 溶解性  |                                                                                                                                                                |
| 纯度   |                                                                                                                                                                |