

# FITC-Dextran (MW 3000-5000)

货号: **B28343**

## FITC-Dextran (MW 3000-5000)

产品信息

生物活性	FITC-Dextran (MW 3000-5000) is a biochemical reagent that can be used as a biological material or organic compound for life science related research.
CAS	60842-46-8
中文名称	
分子量	3000-5000(Average)
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	