

# Phospho-NF-kB p65 (S529) Rabbit mAb

货号: B28415

产品信息

|       |                                                                                               |
|-------|-----------------------------------------------------------------------------------------------|
| 反应    | Human                                                                                         |
| 宿主    | Rabbit                                                                                        |
| 克隆性   | Monoclonal                                                                                    |
| 预测反应  |                                                                                               |
| 应用    | WB FC IP                                                                                      |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>FC:</b> 1:20 - 1:50<br><b>IP:</b> 1:20 - 1:50                 |
| 理论分子量 | 58kDa/59kDa/60kDa                                                                             |
| 实测分子量 | 65kDa                                                                                         |
| 形式    | Liquid                                                                                        |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                  |
| 阳性对照  | HeLa,293T,MCF7,Mouse lung,Mouse kidney,Rat lung,Rat kidney                                    |
| 细胞定位  | Cytoplasm,Nucleus                                                                             |
| 纯化    | Affinity purification                                                                         |

抗原信息


|      |                                                                                                                                         |
|------|-----------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human Phospho-NF-kB p65 (S529).                                                             |
| 序列   | RSTDTTKTHPTIKINGYTGP GTVRISLVTKDPPHRPHPELVGKDCRDGFYEALCPDRCIHSFQNLGIQCVKKRDLEQAIS<br>QRIQTNNNPFQVPIEEQRGDYDLNAVRLCFQVTVRDPSSGRPLRLPPVLS |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | NF-kappa-B is a ubiquitous transcription factor involved in several biological processes. It is held in the cytoplasm in an inactive state by specific inhibitors. Upon degradation of the inhibitor, NF-kappa-B moves to the nucleus and activates transcription of specific genes. NF-kappa-B is composed of NFKB1 or NFKB2 bound to either REL, RELA, or RELB. The most abundant form of NF-kappa-B is NFKB1 complexed with the product of this gene, RELA. Four transcript variants encoding different isoforms have been found for this gene. |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                               |
|-------|-------------------------------|
| 基因ID  | 5970                          |
| 基因名   | RELA                          |
| Swiss | Q04206                        |
| 别名    | NFKB3;p65;NF-kB p65;RELA;CMCU |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)