

NFIB / NF1B2 Rabbit mAb

货号: B28503

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	22kDa/33kDa/47kDa/53kDa/55kDa
实测分子量	47kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human NFIB / NF1B2.
序列	SRTWHERDQDMSSPTTMKKPEKPLFSSASPQDSSPRLSTFPQHHPGIPGVAHSVISTRTPPPSPLPFPTQAILPPAPSSYFSHTIRYPHLPNPQDTLKNYVPSYDPSSPQTSQSWYLG

靶点信息

研究背景	Transcriptional activator of GFAP, essential for proper brain development. Recognizes and binds the palin dromic sequence 5'-TTGGCNNNNNGCCAA-3' present in viral and cellular promoters and in the origin of replication of adenovirus type 2. These proteins are individually capable of activating transcription and replication.
------	---

基因ID	4781
基因名	NFIB
Swiss	O00712
别名	NFIB;CTF;HMGIC/NFIB;NF-I/B;NF1-B;NFI-B;NFI-RED;NFIB2;NFIB3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn