

MiTF Rabbit mAb

货号: B28565

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	40kDa/46kDa/52kDa/55kDa/56kDa/57kDa/58kDa
实测分子量	52kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,K-562,Mouse kidney
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human MiTF.
序列	MLEMLEYNHYQVQTHLENPTKYHIQQAQRQVKQYLSTTLANKHANQVLSLPCPNQPGDHVMPPVPGSSAPNSPMAML TLNSNCEKEGFYKFEEQNRAESECpgmNTHSRASCMQMDDVIDDIISLESSYNEEILGLMDPALQMANTLPVSGNLIDLYG NQGLPPPGLTISNSCPANLPNIKRELTESEARALAKERQKKDNHNLIERRRRFNINDRIKELGTLPKSNPDPMRWNKGTLK ASVDYIRKLQREQQRAKELENRQKKLEHANRHLLRIQ

靶点信息

研究背景	This gene encodes a transcription factor that contains both basic helix-loop-helix and leucine zipper structural features. It regulates the differentiation and development of melanocytes retinal pigment epithelium and is also responsible for pigment cell-specific transcription of the melanogenesis enzyme genes. Heterozygous mutations in the this gene cause auditory-pigmentary syndromes, such as Waardenburg syndrome type 2 and Tietz syndrome. Alternatively spliced transcript variants encoding different isoforms have been identified.
------	---

基因ID	4286
基因名	MITF
Swiss	O75030
别名	MITF;CMM8;COMMAD;MI;WS2;WS2A;bHLHe32

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn