

SERPING1 Rabbit mAb

货号: B28575

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	WB: 1:500 - 1:2000 IP: 1:20 - 1:50
理论分子量	49kDa/55kDa
实测分子量	100kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Human plasma
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human SERPING1.
序列	KDFTCVHQALKGFTTKGVTSVSQIFHSPDLAIRDTFVNASRTLYSSSPRVLSNNSDANLELINTWVAKNTNNKISRLLDSLP SDTRLVLLNAIYLSAKWKTTFDPKKTRMEPFHFKNVSVIKVPMMNSSKKYPVAHFIDQTLKAKVGQLQLSHNLSLVILVPQNLK HRLDMEQALSPSVFKAIMEKLEMSKFQPTLLTPRIKVTTSQDMLSIMEKLEFFDFSVDLNLCLGTEDPDLQVSAMQHQT VLELTETGVEAAAASISVARTLLVFEVQQPFLFVLWDQQHKFPVFMGRVYDPRA

靶点信息

研究背景	This gene encodes a highly glycosylated plasma protein involved in the regulation of the complement cascade. Its protein inhibits activated C1r and C1s of the first complement component and thus regulates complement activation. Deficiency of this protein is associated with hereditary angioneurotic oedema (HANE). Alternative splicing results in multiple transcript variants encoding the same isoform.
------	---

基因ID	710
基因名	SERPING1
Swiss	P05155
别名	SERPING1;C1IN;C1INH;C1NH;HAE1;HAE2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn