

CDX2 Rabbit mAb

货号: B28579

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	33kDa
实测分子量	40kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW620,MCF7
细胞定位	Nucleus
纯化	Affinity purification

抗原信息



抗原信息	Recombinant fusion protein corresponding to Human CDX2.
序列	MYVSYLLDKDVSMYPSSVRHSGGLNLAPQNFVSPQYPDYGGYHVAAAAAAAAANLDSAQSPGPSWPAAYGAPLREDW NGYAPGGAAAAANAVAHGLNGGSPAAAMGYSSPADYHPHHHPHHHPHAAAPSCASGLLQTLNPGPPGAATAAAEQ LSPGGQRRNLCEWMRKPAQQLGSQ

靶点信息

研究背景	This gene is a member of the caudal-related homeobox transcription factor gene family. The encoded protein is a major regulator of intestine-specific genes involved in cell growth and differentiation. This protein also plays a role in early embryonic development of the intestinal tract. Aberrant expression of this gene is associated with intestinal inflammation and tumorigenesis.
------	--

基因ID	1045
基因名	CDX2
Swiss	Q99626
别名	CDX2;CDX-3;CDX2/AS;CDX3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn