

CLASP1 Rabbit mAb

货号: B28648

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	162kDa/163kDa/164kDa/169kDa
实测分子量	169kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,MCF7,Mouse brain,Rat brain
细胞定位	Chromosome,Cytoplasm,Golgi apparatus,centromere,centrosome,cytoskeleton,kinetochore,microtubule o rganizing center,spindle,trans-Golgi network
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human CLASP1.
序列	IMKTLEAHKDSHKEVVRAAEEAASTLASSIHPEQCIKVLCP IIQTADYPINLAAIKMQTKVVERIAKESLLQLLVDIIPG LLQGYD NTESSVRKASVFLVAIYSVIGEDLKPHLAQLTGSKMKLLNLYIKRAQTTNSNSSSSSDVSTHS

靶点信息

研究背景	CLASPs, such as CLASP1, are nonmotor microtubule-associated proteins that interact with CLIPs (e.g., CLI P170; MIM 179838). CLASP1 is involved in the regulation of microtubule dynamics at the kinetochore and throughout the spindle (Maiato et al., 2003 [PubMed 12837247]).[supplied by OMIM, Mar 2008]
------	---

基因ID	23332
基因名	CLASP1
Swiss	Q7Z460
别名	CLASP1;MAST1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn