

DDAH2 Rabbit mAb

货号: B28686

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	WB: 1:500 - 1:2000 IP: 1:20 - 1:50
理论分子量	29kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW620,BT-474,MCF7,Mouse lung,Rat brain
细胞定位	Cytoplasm,Mitochondrion
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human DDAH2.
序列	MGTPGEGLGRCSHALIRGVPESLASGEGAGAGLPALDLAKAQREHGVLGGLRQRLGLQLLELPPEESLPLGPLLGDTAVI QGDTALITRPWSPARRPEVDGVRKALQDLGLRIVEIGDENATLDGTDVLTGREFFVGLSKWTNHRGAEIVADTFRDFAV STVPVSGPSHLRGLCGMGGPRTVVAGSSDAAQKAVRAMAVLTDHPYASLTLPDDAAADCLFLRPGLPGVPPFLLHRGGG DLPNSQEALQKLSDVTLVPVSCSELEKAGAGLSSLCLVLSTRPHS

靶点信息

研究背景	This gene encodes a dimethylarginine dimethylaminohydrolase. The encoded enzyme functions in nitric oxide generation by regulating the cellular concentrations of methylarginines, which in turn inhibit nitric oxide synthase activity. The protein may be localized to the mitochondria. Alternative splicing resulting in multiple transcript variants.
------	--

基因ID	23564
基因名	DDAH2
Swiss	O95865
别名	DDAH2;DDAH;DDAHII;G6a;HEL-S-277;NG30

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn