

# ITCH Rabbit mAb

货号: B28688

## 产品信息

|       |                                                                                                     |
|-------|-----------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                     |
| 宿主    | Rabbit                                                                                              |
| 克隆性   | Monoclonal                                                                                          |
| 预测反应  |                                                                                                     |
| 应用    | WB                                                                                                  |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                           |
| 理论分子量 | 103kDa                                                                                              |
| 实测分子量 | 103kDa                                                                                              |
| 形式    | Liquid                                                                                              |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.       |
| 偶联物   | Unconjugated                                                                                        |
| 阳性对照  | HeLa,293T,Raji,NIH/3T3,C6,Mouse lung                                                                |
| 细胞定位  | cell cortex,cytoplasm,cytoplasmic vesicle,cytosol,extracellular exosome,nucleoplasm,plasma membrane |
| 纯化    | Affinity purification                                                                               |

## 抗原信息


|      |                                                                                                       |
|------|-------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human ITCH .                                              |
| 序列   | MSDSGSQLGSMGSLTMKSQLQITVISAKLKENKKNWFGPSPYVEVTVDGQSKKTEKCNNNTNSPKWKQPLTVIVTPVSKLHFRVWSHQTLKSDVLLGTAAL |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a member of the Nedd4 family of HECT domain E3 ubiquitin ligases. HECT domain E3 ubiquitin ligases transfer ubiquitin from E2 ubiquitin-conjugating enzymes to protein substrates, thus targeting specific proteins for lysosomal degradation. The encoded protein plays a role in multiple cellular processes including erythroid and lymphoid cell differentiation and the regulation of immune responses. Mutations in this gene are a cause of syndromic multisystem autoimmune disease. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Mar 2012] |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                          |
|-------|--------------------------|
| 基因ID  | 83737                    |
| 基因名   | ITCH                     |
| Swiss | Q96J02                   |
| 别名    | ADMFD; AIF4; AIP4; NAPP1 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)