

# FOXA1 Rabbit mAb

货号: B28724

产品信息

|       |                                                                                               |
|-------|-----------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                               |
| 宿主    | Rabbit                                                                                        |
| 克隆性   | Monoclonal                                                                                    |
| 预测反应  |                                                                                               |
| 应用    | WB IHC IF/ICC                                                                                 |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200<br><b>IF/ICC:</b> 1:50 - 1:200          |
| 理论分子量 | 45kDa/49kDa                                                                                   |
| 实测分子量 | 49kDa                                                                                         |
| 形式    | Liquid                                                                                        |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                  |
| 阳性对照  | Rat testis                                                                                    |
| 细胞定位  | Nucleus                                                                                       |
| 纯化    | Affinity purification                                                                         |

抗原信息


|      |                                                                                                                                                                                                                        |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human FOXA1.                                                                                                                                                               |
| 序列   | QKRFKCEKQPGAGGGGGSGGGSGAKGGPESRKDPGASNPADSPLHRGVHGKTGQLEGAPAPGPAASPQTLDHSG<br>ATATGGASELKTPASSTAPPISSGPGALASVPASHPAHGLAPHESQLHLKGDPHYSFNHPFSINNLMSSSEQQHKLDFKAY<br>EQALQYSPYGSTLPASLPLGSASVTRSPIEPSALEPAYYQGVYSRPVLNTS |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                    |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a member of the forkhead class of DNA-binding proteins. These hepatocyte nuclear factors are transcriptional activators for liver-specific transcripts such as albumin and transthyretin, and they also interact with chromatin. Similar family members in mice have roles in the regulation of metabolism and in the differentiation of the pancreas and liver. |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                   |
|-------|-------------------|
| 基因ID  | 3169              |
| 基因名   | FOXA1             |
| Swiss | P55317            |
| 别名    | FOXA1;HNF3A;TCF3A |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)