

DR5 Rabbit mAb

货号: B28758

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP FC
推荐浓度	WB: 1:500 - 1:2000 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	12kDa/45kDa/47kDa
实测分子量	40,48kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human DR5.
序列	RSSQRPGAEDNVLNEIVSILQPTQVPEQEMEVEQPAEPTGVNMLSPGESEHLLPEAEERSQRRRLLVPANEGDPTETLRQ CFDDFADLVPFDSWEPLMRKLGLMDNEIKVAKAEAAAGHRDTLYTMLIKWVNKTGRDASVHTLLDALETLGERLAKQKIED HLLSSGKFMYLEGNADSAMS

靶点信息

研究背景	The protein encoded by this gene is a member of the TNF-receptor superfamily, and contains an intracellular death domain. This receptor can be activated by tumor necrosis factor-related apoptosis inducing ligand and (TNFSF10/TRAIL/APO-2L), and transduces an apoptosis signal. Studies with FADD-deficient mice suggested that FADD, a death domain containing adaptor protein, is required for the apoptosis mediated by this protein. Two transcript variants encoding different isoforms and one non-coding transcript have been found for this gene.
------	---

基因ID	8795
基因名	TNFRSF10B
Swiss	O14763
别名	CD262;DR5;KILLER;KILLER/DR5;TRAIL-R2;TRAILR2;TRICK2;TRICK2A;TRICK2B;TRICKB;ZTNFR9;TNFRSF10B

产品验证



Western blot analysis of DR5 expressed in Jurkat using DR5 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn