

# JAK3 Rabbit mAb

货号: B28766

## 产品信息

|       |                                                                                               |
|-------|-----------------------------------------------------------------------------------------------|
| 反应    | Human                                                                                         |
| 宿主    | Rabbit                                                                                        |
| 克隆性   | Monoclonal                                                                                    |
| 预测反应  |                                                                                               |
| 应用    | WB ICC                                                                                        |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>ICC:</b> 1:50 - 1:200                                         |
| 理论分子量 | 68kDa/121kDa/125kDa                                                                           |
| 实测分子量 | 125kDa                                                                                        |
| 形式    | Liquid                                                                                        |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                  |
| 阳性对照  | Jurkat                                                                                        |
| 细胞定位  | Cytoplasm,Endomembrane system,Peripheral membrane protein                                     |
| 纯化    | Affinity purification                                                                         |

## 抗原信息


|      |                                                                                                                                                                                              |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human JAK3.                                                                                                                                      |
| 序列   | RLTTDSQHFFCKEVAPPRLLEEVAEQCHGPITLDFAINKLKTGGSRPGSYVLRSPQDFDSFLLTVCVQNPLGPDYKGCLIR<br>RSPTGTFLLVGLSRPHSSLRELLATCWDGGLHVDGVAVTLTSCCIPRPKEKSNLIVVQRGHSPPTSSLVQPQSQYQLSQMT<br>FHKIPADSLEWHENLGHG |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | The protein encoded by this gene is a member of the Janus kinase (JAK) family of tyrosine kinases involved in cytokine receptor-mediated intracellular signal transduction. It is predominantly expressed in immune cells and transduces a signal in response to its activation via tyrosine phosphorylation by interleukin receptors. Mutations in this gene are associated with autosomal SCID (severe combined immunodeficiency disease). |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                                       |
|-------|---------------------------------------|
| 基因ID  | 3718                                  |
| 基因名   | JAK3                                  |
| Swiss | P52333                                |
| 别名    | JAK3;JAK-3;JAK3_HUMAN;JAKL;L-JAK;LJAK |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)