

Asparagine synthetase Rabbit mAb

货号: B28802

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	54kDa/62kDa/64kDa
实测分子量	61kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SGC-7901,U-87MG,Mouse brain,Mouse pancreas,Mouse spleen,Rat testis
细胞定位	cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Asparagine synthetase.
序列	DSSLVAATLLKQLKEAQVQYPLQTFAIGMEDSPDLLAARKVADHIGSEHYEVLFNSEEGIQALDEVIFSLETYDITTVRASVG MYLISKYIRKNTDSVVIFSGEGSDELTQGYIYFHKAPSPEKAEESERLLRELYLFDVLRADRTTAAHGLELRVPFLDHRFSSY YLSLPPEMRIPKNGIEKHLLRETFEDSNLIPKEILWRPKEAFSDGITSVKNSWFKILQEYVEHQVDDAMMANAAQKFPFNT P KTKEGYYYRQVFERHYPGRADWLSHYWMPKWINATDPSARTLTHYKSAVKA

靶点信息

研究背景	The protein encoded by this gene is involved in the synthesis of asparagine. This gene complements a mu tation in the temperature-sensitive hamster mutant ts11, which blocks progression through the G1 phase of the cell cycle at nonpermissive temperature. Alternatively spliced transcript variants have been describ ed for this gene.
------	---

基因ID	440
基因名	ASNS
Swiss	P08243
别名	ASNS;ASNSD;TS11

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn