

CD79a Rabbit mAb

货号: B28811

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	45-55kDa
实测分子量	45-55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji,Daudi,Mouse spleen,Rat spleen,Rat thymus
细胞定位	Cell membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human CD79a.
序列	FRKRWQNEKLGLDAGDEYEDENLYEGLNLDDCSMYEDISRGLQGTYQDVGSLNIGDVQLEK

靶点信息

研究背景	The B lymphocyte antigen receptor is a multimeric complex that includes the antigen-specific component , surface immunoglobulin (Ig). Surface Ig non-covalently associates with two other proteins, Ig-alpha and I g-beta, which are necessary for expression and function of the B-cell antigen receptor. This gene encodes the Ig-alpha protein of the B-cell antigen component. Alternatively spliced transcript variants encoding dif ferent isoforms have been described. [provided by RefSeq, Jul 2008]
------	---

基因ID	973
基因名	CD79A
Swiss	P11912
别名	CD79A;IGA;MB-1;CD79a molecule;CD79a

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn