

P70 S6 Kinase beta Rabbit mAb

货号: B28838

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	16kDa/53kDa
实测分子量	53kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,SH-SY5Y,HL-60,SW620,MCF7,Mouse liver,Mouse kidney,Mouse spleen,Rat liver
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human P70 S6 Kinase beta.
序列	MAAVFDLDLETEEGSEGEPELSPADACPLAELRAAGLEPVGHYEELTETSVNVGPERIGPHCFELLRVLGKGGYGKV FQVRKVQGNTNLGKIYAMKVLRKAKIVRNAKDTAHTRAERNILESVKHPFIVELAYAFQTGGKLYLILECLSGGELFTHLEREGI FLEDT

靶点信息

研究背景	This gene encodes a member of the RSK (ribosomal S6 kinase) family of serine/threonine kinases. This kinase contains a kinase catalytic domain and phosphorylates the S6 ribosomal protein and eukaryotic translation initiation factor 4B (eIF4B). Phosphorylation of S6 leads to an increase in protein synthesis and cell proliferation.
------	---

基因ID	6199
基因名	RPS6KB2
Swiss	Q9UBS0
别名	RPS6KB2;KLS;P70-beta;P70-beta-1;P70-beta-2;S6K-beta2;S6K2;SRK;STK14B;p70(S6K)-beta

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn