

IFNAR1 Rabbit mAb

货号: B28858

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP FC
推荐浓度	WB: 1:500 - 1:2000 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	49kDa/56kDa/63kDa
实测分子量	90,130kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SP2/0,HeLa,A-549
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human IFNAR1.
序列	PSLKPSSSIDEYFSEQPLKNLLSTSEEQIEKCFIENISTIATVEETNQTDDEHKKYSSQTSQDSGNYSNEDESESKTSEELQ QDFV

靶点信息

研究背景	The protein encoded by this gene is a type I membrane protein that forms one of the two chains of a receptor for interferons alpha and beta. Binding and activation of the receptor stimulates Janus protein kinases, which in turn phosphorylate several proteins, including STAT1 and STAT2. The encoded protein also functions as an antiviral factor.
基因ID	3454

基因名	IFNAR1
Swiss	P17181
别名	IFNAR1;AVP;IFN-alpha-REC;IFNAR;IFNBR;IFRC

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn