

# FOXO3A Rabbit mAb

货号: B28885

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	48kDa/71kDa
实测分子量	90kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Cytoplasm,Nucleus,cytosol
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein corresponding to Human FOXO3A.
序列	AVSMD

## 靶点信息

研究背景	This gene belongs to the forkhead family of transcription factors which are characterized by a distinct fork head domain. This gene likely functions as a trigger for apoptosis through expression of genes necessary for cell death. Translocation of this gene with the MLL gene is associated with secondary acute leukemia. Alternatively spliced transcript variants encoding the same protein have been observed.
基因ID	2309

基因名	FOXO3
Swiss	O43524
别名	FOXO3;AF6q21;FKHRL1;FKHRL1P2;FOXO2;FOXO3A;FoxO3a

产品验证



Western blot analysis of FOXO3A expressed in 293T using FOXO3A Rabbit mAb at 1:1000. Secondary anti body: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% no n-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)