

ROCK2 Rabbit mAb

货号: B28887

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	160kDa
实测分子量	161kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,THP-1, Mouse lung,Mouse brain,Rat brain
细胞定位	Cell membrane,Cytoplasm,Nucleus,Peripheral membrane protein,centrosome,cytoskeleton,microtubule o rganizing center
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human ROCK2 .
序列	LDSKDSIEQLRSQLQALHIGLDSSSIGSGPGDAEADDGFPESRLEGWLSLPVRNNTKKFGWVKKYVIVSSKKILFYDSEQ DKEQSNPYMVLIDIDKLFHVRPVTQTDVYRADAKEIPRIFQILYANEGESKKEQEFPPVPGVGEKSNIYCHKGHEFIPTLYHFPT NCEACMKPLWHMFKPPPALECCRCHIKCHKDHMDKKEEIIAPCKVYYDISTAKNLLLLANSTEEQQKWVSRLVKKIPKPPA PDPFARSSPRTSMKIQQNQSIRRPQRQLAPNKPS

靶点信息

研究背景	The protein encoded by this gene is a serine/threonine kinase that regulates cytokinesis, smooth muscle contraction, the formation of actin stress fibers and focal adhesions, and the activation of the c-fos serum response element. This protein, which is an isozyme of ROCK1 is a target for the small GTPase Rho.
------	---

基因ID	9475
基因名	ROCK2
Swiss	O75116
别名	ROCK2;ROCK-II

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn