

AHCY Rabbit mAb

货号: B28890

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	44kDa/47kDa
实测分子量	47kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Mouse pancreas,Rat liver
细胞定位	Cytoplasm,Melanosome
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human AHCY.
序列	GILKVPAINVNDsvTKSKFDnLYGCRSLIDGIKRATDVMiAGKVAVVAGYGDVGKGCAQALRGFGARVIITEIDPINALQAA MEGYEVTTmDEACqEGNIFVTTTGCIDIILGRHFEQMKDDAIVCNIGHFDVEIDVKWLNENAVEKVNIKPQVDryRLKNGR RIILLAEGRlVNLGCAMGHPSFVMSNSFTnQVMAQIELWTHPDkYPVGvHFLPKKLDEAVAEaHLGKLNVKLTKLteKQAQ YLGmSCDGPfKPDHYR

靶点信息

研究背景	S-adenosylhomocysteine hydrolase belongs to the adenosylhomocysteinase family. It catalyzes the reversible hydrolysis of S-adenosylhomocysteine (AdoHcy) to adenosine (Ado) and L-homocysteine (Hcy). Thus , it regulates the intracellular S-adenosylhomocysteine (SAH) concentration thought to be important for transmethylation reactions. Deficiency in this protein is one of the different causes of hypermethioninemia. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	191
基因名	AHCY
Swiss	P23526
别名	AHCY;SAHH;adoHcyase

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn