

GPRC5B Rabbit mAb

货号: B28917

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	52kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,BxPC-3,K-562,Mouse brain,Mouse kidney,Rat brain
细胞定位	cell surface,cytosol,extracellular exosome,extracellular space,nucleolus,nucleoplasm,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human GPRC5B.
序列	TPNYFDTSQPRMRETAFEEDVQLPRAYMENKAFSMDEHNAAALRTAGFPNGSLGKRPSGSLGKRPSAPFRSNVYQPTEMA VVLNGGTIPTAPPSHTGRHLW

靶点信息

研究背景	This gene encodes a member of the type 3 G protein-coupled receptor family. Members of this superfamil y are characterized by a signature 7-transmembrane domain motif. The encoded protein may modulate i nsulin secretion and increased protein expression is associated with type 2 diabetes. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2015]
基因ID	51704

基因名	GPRC5B
Swiss	Q9NZH0
别名	

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn