

# LSP1 Rabbit mAb

货号: B28936

产品信息

|       |                                                                                               |
|-------|-----------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                               |
| 宿主    | Rabbit                                                                                        |
| 克隆性   | Monoclonal                                                                                    |
| 预测反应  |                                                                                               |
| 应用    | WB IHC IF/ICC                                                                                 |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200<br><b>IF/ICC:</b> 1:50 - 1:200          |
| 理论分子量 | 30kDa/37kDa/50kDa                                                                             |
| 实测分子量 | 52kDa                                                                                         |
| 形式    | Liquid                                                                                        |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                  |
| 阳性对照  | Mouse lung,Mouse spleen                                                                       |
| 细胞定位  | Cell membrane,Cytoplasmic side,Peripheral membrane protein                                    |
| 纯化    | Affinity purification                                                                         |

抗原信息


|      |                                                                                                                                                                                                                                                                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human LSP1.                                                                                                                                                                                                                                                                                                              |
| 序列   | MAEASSDPGAEEEREELLGPTAQWSVEDEEEEAVHEQCQHERDRQLQAQDEEGGGHVPERPKQEMLLSLKPSEAPELDED<br>EGFGDWSQRPEQRQQHEGAQGALDSGEPPQCRSPEGEQEDRPGLHAYEKEDSDEVHLEELSLSKEGPGPEDTVQDNL<br>GAAGAEEEQEEHQKCQQRTPSPPLVLEGTIEQSSPPLSPTTKLIDRTESLNRSIEKSNSVKKSQPDLPISKIDQWLEQYTQAI<br>ETAGRTPKLARQASIELPSMAVASTKSRWETGEVQAQSAAKTPSCKDIVAGDMSKSLWEQKGGSKTSSTIKSTPSGKRY<br>KFVATGHGKYEKVLVEGGPAP |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                        |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes an intracellular F-actin binding protein. The protein is expressed in lymphocytes, neutrophils, macrophages, and endothelium and may regulate neutrophil motility, adhesion to fibrinogen matrix proteins, and transendothelial migration. Alternative splicing results in multiple transcript variants encoding different isoforms. |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                |
|-------|----------------|
| 基因ID  | 4046           |
| 基因名   | LSP1           |
| Swiss | P33241         |
| 别名    | LSP1;WP34;pp52 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)