

ARRB1 Rabbit mAb

货号: **B28949**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	46kDa/47kDa
实测分子量	50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,A-549,Mouse lung,Rat lung,Rat brain
细胞定位	Cell membrane,Cell projection,Cytoplasm,Cytoplasmic vesicle,Membrane,Nucleus,clathrin-coated pit,pseudopodium
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human ARRB1.
序列	RKVQYAPERPGPQPTAETTRQFLMSDKPLHLEASLDKEIYYHGEPISVNVHVTNNTNKT VKKIKISVRQYADICLFNTAQYKCPVAMEEADDTVAPSSTFCVKVYTLTPFLANNREKRG LALDGLKHEDTNLASSTLLREGANREILGIIVSYKVKV KLVVSRGGLLGLASSDVAVELPFTLMHPKPKEEPPHREVPENETPVDTNLIELDTNDDDIVFEDFARQRLKGMKDDKEEEDGTGSPQLNNR

靶点信息

研究背景	Members of arrestin/beta-arrestin protein family are thought to participate in agonist-mediated desensitization of G-protein-coupled receptors and cause specific dampening of cellular responses to stimuli such as hormones, neurotransmitters, or sensory signals. Arrestin beta 1 is a cytosolic protein and acts as a cofactor in the beta-adrenergic receptor kinase (BARK) mediated desensitization of beta-adrenergic receptors. Besides the central nervous system, it is expressed at high levels in peripheral blood leukocytes, and thus the BARK/beta-arrestin system is believed to play a major role in regulating receptor-mediated immune functions. Alternatively spliced transcripts encoding different isoforms of arrestin beta 1 have been described.
基因ID	408
基因名	ARRB1
Swiss	P49407
别名	ARRB1;ARB1;ARR1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn