

PEDF Rabbit mAb

货号: B28950

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC
推荐浓度	WB: 1:500 - 1:2000 FC: 1:20 - 1:50
理论分子量	46kDa
实测分子量	46kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Human plasma,Mouse liver,Rat liver
细胞定位	Melanosome,Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human PEDF.
序列	PDEISILLLGVAHFKGQWVTKFDSRKTSLEDFYLDEERTVRVPMMSDPKAVLRYGLDSDLCKIAQLPLTGSMsIIFFLPLKV TQNLTLIEESLTSEFIHD

靶点信息

研究背景	The protein encoded by this gene is a member of the serpin family, although it does not display the serin e protease inhibitory activity shown by many of the other serpin family members. The encoded protein is secreted and strongly inhibits angiogenesis. In addition, this protein is a neurotrophic factor involved in n eural differentiation in retinoblastoma cells.
基因ID	5176

基因名	SERPINF1
Swiss	P36955
别名	SERPINF1;EPC-1;OI12;OI6;PEDF;PIG35

产品验证



Western blot analysis of PEDF expressed in Human plasma,Mouse liver,Rat liver using PEDF Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn