

GLUR Rabbit mAb

货号: B28957

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC
推荐浓度	WB: 1:500 - 1:2000 FC: 1:20 - 1:50
理论分子量	47kDa/50kDa/51kDa/53kDa/56kDa
实测分子量	56kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,HeLa,A-549,Jurkat,Mouse brain,Mouse heart,Mouse kidney
细胞定位	Cytoplasm,Mitochondrion
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human GLUR.
序列	MISTNCTEELENAGVEVLKFSQVKEVKKTLSGLEVSMVTAVPGRLPVMTMIPDVDCLLWAIGRVPNTKDLSLNKLGITDD KGHIIVDEFQNTNVKGIYAVGDVCGKALLTPVAIAAGRKLAHRLFYKEDSKLDYNNIPTVVFSHPPIGTVGLTEDEAIHKYGI ENVKTYSTSFTPMYHAVTKRKTKCVMKMCANKEEKVVGIIHQGLGCDEMLQGFVAVKMGATKADFDNTVAIHPTSSE ELVTLR

靶点信息

研究背景	This gene encodes a member of the class-I pyridine nucleotide-disulfide oxidoreductase family. This enzyme is a homodimeric flavoprotein. It is a central enzyme of cellular antioxidant defense, and reduces oxidized glutathione disulfide (GSSG) to the sulfhydryl form GSH, which is an important cellular antioxidant. Rare mutations in this gene result in hereditary glutathione reductase deficiency. Multiple alternatively spliced transcript variants encoding different isoforms have been found.
------	--

基因ID	2936
基因名	GSR
Swiss	P00390
别名	GSR;HEL-75;HEL-S-122m

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn