

Prolactin Rabbit mAb

货号: B28961

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	25kDa
实测分子量	23kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Prolactin.
序列	LPICPGGAARCQVTLRDLFDRAVVLSHYIHNLSEMSEFDKRYTHGRGFITKAINSCHTSSLATPEDKEQAQQMNQKDFL SLIVSILRSWNEPLYHLVTEVRGMQEAPAILSKAVEIEEQTKRLLEGMELIVSQVHPETKENEIYPVWSGLPSLQMADEESRL SAYYNLLHCLRRDSHKIDNYLKLLKCRIIHNNNC

靶点信息

研究背景	This gene encodes the anterior pituitary hormone prolactin. This secreted hormone is a growth regulator for many tissues, including cells of the immune system. It may also play a role in cell survival by suppressing apoptosis, and it is essential for lactation. Alternative splicing results in multiple transcript variants that encode the same protein.
基因ID	5617

基因名	PRL
Swiss	P01236
别名	PRL;GHA1;prolactin

产品验证



Western blot analysis of Prolactin expressed in Mouse brain,Rat brain using Prolactin Rabbit mAb at 1:100 0. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Block ing buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn