

PLTP Rabbit mAb

货号: B29015

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	44kDa/49kDa/54kDa
实测分子量	80kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,SW620
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human PLTP.
序列	VLLNSLLDTVPVRSSVDELVGIDYSLMKDPVASTSNLDMDFRGAFFPLTERNWSLPNRAVEPQLQEEERMVYVAFSEFFF DSAMESYFRAGALQLLLVGDKVPHDLDMLLRATYFGSIVLLSPAVIDSPLKLELRVLAPPRCTIKPSGTTISVTASVTIALVPP DQPEVQLSSMTMDARLSAKMALRGKALRTQLDLRRFRIYSNHSALSLALIPLQAPLKTMLQIGVMPMLNERTWRGVQIPL PEGINFVHEVVTNHAGFLTIGADLHFAKGLREVIEKNRPADVRASTAPTSTAAS

靶点信息

研究背景	The protein encoded by this gene is one of at least two lipid transfer proteins found in human plasma. Th e encoded protein transfers phospholipids from triglyceride-rich lipoproteins to high density lipoprotein (H DL). In addition to regulating the size of HDL particles, this protein may be involved in cholesterol metabo lism. At least two transcript variants encoding different isoforms have been found for this gene.
基因ID	5360

基因名	PLTP
Swiss	P55058
别名	PLTP;BPIFE;HDLCQ9

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn