

Cytochrome C Rabbit mAb

货号: B29017

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	11kDa
实测分子量	15kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,SH-SY5Y,BT-474
细胞定位	Mitochondrion intermembrane space
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Cytochrome C.
序列	MGDVEKGKKIFIMKCSQCHTVEKGGKHKTGPNLHGLFGRKTGQAPGYSYAANKNKGIIWGEDTLMEYLENPKKYIPGTK MIFVGIKKKEERADLIAYLK KATNE

靶点信息

研究背景	This gene encodes a small heme protein that functions as a central component of the electron transport c hain in mitochondria. The encoded protein associates with the inner membrane of the mitochondrion whe re it accepts electrons from cytochrome b and transfers them to the cytochrome oxidase complex. This pr otein is also involved in initiation of apoptosis. Mutations in this gene are associated with autosomal domi nant nonsyndromic thrombocytopenia. Numerous processed pseudogenes of this gene are found through out the human genome.
------	---

基因ID	54205
基因名	CYCS
Swiss	P99999
别名	CYCS;CYC;HCS;THC4;cytochrome c

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn