

CD40L Rabbit mAb

货号: B29113

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	29kDa
实测分子量	36kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Molt-4,Mouse thymus
细胞定位	Cell membrane,Secreted,Single-pass type II membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human CD40L.
序列	MQKGDQNPQIAAHVISEASSKTTSVLQWAEKGYTMSNNLVTLENGKQLTVKRQGLYYIAQVTFCSNREASSQAPFIASL CLKSPGRFERILLRAANTHSSAKPCGQQSIHLGGVFELQPGASVFNVTDP SQVSHGTGFTSFGLLKL

靶点信息

研究背景	The protein encoded by this gene is expressed on the surface of T cells. It regulates B cell function by eng aging CD40 on the B cell surface. A defect in this gene results in an inability to undergo immunoglobulin c lass switch and is associated with hyper-IgM syndrome.
基因ID	959
基因名	CD40LG

Swiss	P29965
别名	CD154;CD40L;HIGM1;IGM;IMD3;T-BAM;TNFSF5;TRAP;gp39;hCD40L;CD40LG

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn