

RBP4 Rabbit mAb

货号: B29129

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	23kDa
实测分子量	23kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW620,HeLa,OVCAR3,Mouse liver,Mouse spleen
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human RBP4.
序列	ERDCRVSSFRVKENFDKARFSGTWYAMAKKDPEGLFLQDNIVAEFVSVDETGQMSATAKGRVRLNNWDVCADMVGTF TDTEDPAKFKMKYWGVASFLLQKGNDHWIVDTDYDYAVQYSCRLNLDGTCADSYSFVFSRDPNGLPPEAQKIVRQR QEELCLARQYRLIVHNGYCDGRSERNLL

靶点信息

研究背景	This protein belongs to the lipocalin family and is the specific carrier for retinol (vitamin A alcohol) in the blood. It delivers retinol from the liver stores to the peripheral tissues. In plasma, the RBP-retinol complex interacts with transthyretin which prevents its loss by filtration through the kidney glomeruli. A deficiency of vitamin A blocks secretion of the binding protein posttranslationally and results in defective delivery and supply to the epidermal cells.
------	--

基因ID	5950
基因名	RBP4
Swiss	P02753
别名	RBP4;MCOPCB10;RDCCAS

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn