

APC Rabbit mAb

货号: B29156

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	300kDa/311kDa
实测分子量	160kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW480,Mouse brain,Rat brain
细胞定位	Cell junction,Cell membrane,Cell projection,Cytoplasm,adherens junction,cytoskeleton,lamellipodium,ruffle membrane
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human APC.
序列	QSKENQVSAKGTWRKIKENEFSPNTSTSQTVSSGATNGAESKTLIYQMAPAVSKTEDVWVRIEDCPINNPRSGRSPTGNT PPVIDSVSEKANPNIKDSKDNQAKQNVGNGSVPMRTVGLENRLNSFIQVDAPDQKGTEIKPGQNNPVPVSETNESSIVER T

靶点信息

研究背景	This gene encodes a tumor suppressor protein that acts as an antagonist of the Wnt signaling pathway. It is also involved in other processes including cell migration and adhesion, transcriptional activation, and a poptosis. Defects in this gene cause familial adenomatous polyposis (FAP), an autosomal dominant pre-m alignant disease that usually progresses to malignancy. Disease-associated mutations tend to be clustere d in a small region designated the mutation cluster region (MCR) and result in a truncated protein product .
基因ID	324
基因名	APC
Swiss	P25054
别名	APC;BTPS2;DP2;DP2.5;DP3;GS;PPP1R46

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn