

Apolipoprotein A II Rabbit mAb

货号: B29162

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	11kDa
实测分子量	11kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Human plasma
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Apolipoprotein A II.
序列	MKLLAATVLLLTICSLEGALVRRQAKEPCVESLSQYFQTVTDYGKDLMEKVKSPELQAEAKSYFEKSKEQLTPLIKKAGTEL VNFLSYFVELGTQPATQ

靶点信息

研究背景	This gene encodes apolipoprotein (apo-) A-II, which is the second most abundant protein of the high density lipoprotein particles. The protein is found in plasma as a monomer, homodimer, or heterodimer with a polipoprotein D. Defects in this gene may result in apolipoprotein A-II deficiency or hypercholesterolemia.
基因ID	336

基因名	APOA2
Swiss	P02652
别名	APOA2;Apo-AII;ApoA-II;apoAII

产品验证



Western blot analysis of Apolipoprotein A II expressed in Human plasma using Apolipoprotein A II Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn