

OGG1 Rabbit mAb

货号: B29164

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	22kDa/36kDa/38kDa/39kDa/40kDa/45kDa/47kDa
实测分子量	39kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2
细胞定位	Mitochondrion,Nucleus,Nucleus matrix,Nucleus speckle,nucleoplasm
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human OGG1.
序列	MPARALLPRRMGHRTLASTPALWASIPCRSELRLDLVLPSGQSFRWREQSPAHWSGVLADQVWTLTQTTEEQLHCTVYR GDKSQASRPTPDELEAVRKYFQLDVTLAQLYHHWGSVDSHFQEVAQKFQGVRLLRQDPIECLFSFICSSNNNIARITGMV ERLCQAFGPRLIQLDDVTYHGFPSLQALAGPEVEAHLRKLGLGYRARYVSASARAILEEQGGLAWLQQLRESSYEEAHKAL CILPGVGTKVADCICLMALDKPQAVPVDVHMWHIAQRDYSWHPTTSQAKGSPSQTNKELGNFFRSLWGPYAGWAQAVL FSADLRQSRHAQEPPAKRRKSGSGPEG

靶点信息

研究背景	This gene encodes the enzyme responsible for the excision of 8-oxoguanine, a mutagenic base byproduct which occurs as a result of exposure to reactive oxygen. The action of this enzyme includes lyase activity for chain cleavage. Alternative splicing of the C-terminal region of this gene classifies splice variants into two major groups, type 1 and type 2, depending on the last exon of the sequence. Type 1 alternative splice variants end with exon 7 and type 2 end with exon 8. All variants share the N-terminal region in common, which contains a mitochondrial targeting signal that is essential for mitochondrial localization. Many alternative splice variants for this gene have been described, but the full-length nature for every variant has not been determined.
基因ID	4968
基因名	OGG1
Swiss	O15527
别名	OGG1;HMMH;HOGG1;MUTM;OGH1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn