

NKX2.2 Rabbit mAb

货号: B29173

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	IHC IF/ICC
推荐浓度	IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	30kDa
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	mouse brain,C6,rat brain
细胞定位	nucleoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human NKX2.2.
序列	MSLTNTKTGFSVKDILDLPDTNDEEGSVAEGPEEENEGPEPAKRAGPLGQGALDAVQSLPLKNPFYDSSDNPYTRWLAST EGLQYSLHGLAAGAPPQDSSSKSPEPSADESPDNDKETPGGGGDAGKKRK

靶点信息

研究背景	The protein encoded by this gene contains a homeobox domain and may be involved in the morphogenes is of the central nervous system. This gene is found on chromosome 20 near NKX2-4, and these two gene s appear to be duplicated on chromosome 14 in the form of TITF1 and NKX2-8. The encoded protein is lik ely to be a nuclear transcription factor.
基因ID	4821

基因名	NKX2-2
Swiss	O95096
别名	NKX2-2;NKX2.2;NKX2B

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn