

PIAS1 Rabbit mAb

货号: B29209

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	36kDa/71kDa/72kDa
实测分子量	76kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2 cells, HEK-293 cells, HeLa cells, Jurkat cells
细胞定位	Nucleus,Nucleus speckle,PML body
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human PIAS1.
序列	VQEVASASYNGVDGCLSTLEHQVASHHQSSNKNKKVEVIDLTIDSSSDEEEEEPSAKRTCPSLSPTSPLNNKGILSLPHQA SPVSRTPSLPAVDTSYINTSLIQDYRHFPHMTMPYDLQGLDFFPFLSGDNQHYNTSLAAAAAAAAVSDDQDLLHSSRFFPY TSSQMFLDQLSAGGSTSLPTTNGSSSGSNSSLVSSNSLRSHSHTVTNRSSTDASIFGIIPDIISLD

靶点信息

研究背景	This gene encodes a member of the protein inhibitor of activated STAT (PIAS) family. PIAS proteins function as SUMO E3 ligases and play important roles in many cellular processes by mediating the sumoylation of target proteins. This protein plays a central role as a transcriptional coregulator of numerous cellular pathways including the STAT1 and nuclear factor kappaB pathways. Alternate splicing results in multiple transcript variants.
------	--

基因ID	8554
基因名	PIAS1
Swiss	O75925
别名	PIAS1;DDXBP1;GBP;GU/RH-II;ZMIZ3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn