

AP2 alpha Rabbit mAb

货号: B29217

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	40kDa/47kDa/48kDa
实测分子量	48kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human AP2 alpha.
序列	GLDPRRDYRRHEDLLHGPHALSSGLGDSIHSPLPHAIIEVPHVEDPGINIPDQTVIKKGPVSLSKSNSNAVSAIPINKDNL

靶点信息

研究背景	The protein encoded by this gene is a transcription factor that binds the consensus sequence 5'-GCCNNN GGC-3'. The encoded protein functions as either a homodimer or as a heterodimer with similar family members. This protein activates the transcription of some genes while inhibiting the transcription of others. Defects in this gene are a cause of branchiooculofacial syndrome (BOFS). Three transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	7020
基因名	TFAP2A
Swiss	P05549
别名	TFAP2A;AP-2;AP-2alpha;AP2TF;BOFS;TFAP2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn