

VAMP8 Rabbit mAb

货号: B29219

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	11kDa
实测分子量	15kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,NIH/3T3,PC-12,Mouse kidney,Mouse liver,Rat lung
细胞定位	Cell membrane,Early endosome membrane,Late endosome membrane,Lysosome membrane,Single-pass type IV membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human VAMP8.
序列	MEEASEGGGNDVRNLQSEVEGVKNIMTQNVERILARGENLEHLRNKTEDLEATSEHFKTTSQKVARKFWWKNVK

靶点信息

研究背景	This gene encodes an integral membrane protein that belongs to the synaptobrevin/vesicle-associated m embrane protein subfamily of soluble N-ethylmaleimide-sensitive factor attachment protein receptors (SN AREs). The encoded protein is involved in the fusion of synaptic vesicles with the presynaptic membrane.
------	--

基因ID	8673
基因名	VAMP8
Swiss	Q9BV40
别名	VAMP8;EDB;VAMP-8

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn