

ADAM10 Rabbit mAb

货号: B29248

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	WB: 1:500 - 1:2000 IP: 1:20 - 1:50
理论分子量	49kDa/84kDa
实测分子量	84kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	RAW264.7
细胞定位	Cell membrane,Endomembrane system,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human ADAM10.
序列	TTSAEKNTCQLYIQTDHLFFKYGTREAVIAQISSHVKAIDTIYQTDFSGIRNISFMVKRIRINTTADEKDPTNPFRFPNIGVE KFLELNSEQNHDDYCLAYVFTDRDFDDGVLGLAWVVGAPSGSSGGICEKSKLYSDGKKKSLNTGIITVQNYGSHVPPKVSH ITFAHEVGHNFGSPHDSGTECTPGESKNLGQKENGNYIMYARATSGDKLNNNKFSLCSIRNISQVLEKKRNNCFVESGQPI CGNGMVEQGEECDGYSQCKDECCFDANQPEGRKCKLKPG

靶点信息

研究背景	Members of the ADAM family are cell surface proteins with a unique structure possessing both potential a dhesion and protease domains. This gene encodes and ADAM family member that cleaves many proteins including TNF-alpha and E-cadherin. Alternate splicing results in multiple transcript variants encoding diff erent proteins that may undergo similar processing.
------	---

基因ID	102
基因名	ADAM10
Swiss	O14672
别名	ADAM10;AD10;AD18;CD156c;CDw156;HsT18717;MADM;RAK;kuz

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn